

## Biomimetic materials to characterize bacteria-host interactions

Stones, Daniel; Al-Saedi, Fitua; Vaz, Diana; Perez Soto, Nicolas; Krachler, Anne Marie

DOI:  
[10.3791/53400](https://doi.org/10.3791/53400)

License:  
None: All rights reserved

Document Version  
Peer reviewed version

Citation for published version (Harvard):  
Stones, D, Al-Saedi, F, Vaz, D, Perez Soto, N & Krachler, AM 2015, 'Biomimetic materials to characterize bacteria-host interactions', *Journal of Visualized Experiments*, no. 105, e53400. <https://doi.org/10.3791/53400>

[Link to publication on Research at Birmingham portal](#)

### General rights

Unless a licence is specified above, all rights (including copyright and moral rights) in this document are retained by the authors and/or the copyright holders. The express permission of the copyright holder must be obtained for any use of this material other than for purposes permitted by law.

- Users may freely distribute the URL that is used to identify this publication.
- Users may download and/or print one copy of the publication from the University of Birmingham research portal for the purpose of private study or non-commercial research.
- User may use extracts from the document in line with the concept of 'fair dealing' under the Copyright, Designs and Patents Act 1988 (?)
- Users may not further distribute the material nor use it for the purposes of commercial gain.

Where a licence is displayed above, please note the terms and conditions of the licence govern your use of this document.

When citing, please reference the published version.

### Take down policy

While the University of Birmingham exercises care and attention in making items available there are rare occasions when an item has been uploaded in error or has been deemed to be commercially or otherwise sensitive.

If you believe that this is the case for this document, please contact [UBIRA@lists.bham.ac.uk](mailto:UBIRA@lists.bham.ac.uk) providing details and we will remove access to the work immediately and investigate.

**Biomimetic materials to characterize bacteria-host interactions**

**AUTHORS:**

Stones, Daniel H<sup>\*</sup>  
Institute of Microbiology and Infection  
School of Biosciences  
University of Birmingham  
Birmingham, UK  
[d.stones@bham.ac.uk](mailto:d.stones@bham.ac.uk)

Al-Saedi, Fitua<sup>\*</sup>  
Institute of Microbiology and Infection  
School of Biosciences  
University of Birmingham  
Birmingham, UK  
[fma332@bham.ac.uk](mailto:fma332@bham.ac.uk)

Vaz, Diana  
Institute of Microbiology and Infection  
School of Biosciences  
University of Birmingham  
Birmingham, UK  
[d.p.vaz@bham.ac.uk](mailto:d.p.vaz@bham.ac.uk)

Perez-Soto, Nicolas  
Institute of Microbiology and Infection  
School of Biosciences  
University of Birmingham  
Birmingham, UK  
[nap375@bham.ac.uk](mailto:nap375@bham.ac.uk)

<sup>\*</sup>These authors have contributed equally.

**CORRESPONDING AUTHOR:**

Krachler, Anne M  
Institute of Microbiology and Infection  
School of Biosciences  
University of Birmingham  
Birmingham, UK  
[a.krachler@bham.ac.uk](mailto:a.krachler@bham.ac.uk)

**KEYWORDS:**

Host-pathogen interaction, bacterial adhesion, host cell attachment, adhesin, biomimicry, bioengineering, chemical biology

**SHORT ABSTRACT:**

Bacterial attachment to host cells is a key step during host colonization and infection. This protocol describes the generation of polymer-coupled recombinant adhesins as biomimetic materials which allow analysis of the contribution of individual adhesins to these processes, independent of other bacterial factors.

**LONG ABSTRACT:**

Bacterial attachment to host cells is one of the earliest events during bacterial colonization of host tissues and thus a key step during infection. The biochemical and functional characterization of adhesins mediating these initial bacteria-host interactions is often compromised by the presence of other bacterial factors, such as cell wall components or secreted molecules, which interfere with the analysis. This protocol describes the production and use of biomimetic materials, consisting of pure recombinant adhesins chemically coupled to commercially available, functionalized polystyrene beads, which have been used successfully to dissect the biochemical and functional interactions between individual bacterial adhesins and host cell receptors. Protocols for different coupling chemistries, allowing directional immobilization of recombinant adhesins on polymer scaffolds, and for assessment of the coupling efficiency of the resulting “bacteriomimetic” materials are also discussed. We further describe how these materials can be used as a tool to inhibit pathogen mediated cytotoxicity and discuss scope, limitations and further applications of this approach in studying bacterial - host interactions.

**INTRODUCTION:**

Dissecting the interactions between bacterial adhesins and host surface receptors at the host-pathogen interface is an essential step towards our understanding of the underlying mechanisms driving bacterial adhesion. Ultimately, this will help us to identify strategies to interfere with these processes during infections. In conducting such studies, we often face a dilemma: Biochemical and biophysical analysis of the molecular mechanisms of adhesin binding requires their separation from other cell wall components, which may interfere with adhesion events. On the other hand, the use of recombinant soluble proteins over-simplifies the adhesion event, by disregarding protein anchoring in the bacterial cell wall and multivalency of binding achieved through bacterial surface display. Equally, from the host cell’s perspective, encountering and binding to single adhesin molecules does not always have the same impact in terms of plasma membrane organization, membrane fluidity and receptor clustering<sup>1</sup> and this, ultimately, makes it difficult to evaluate the impact of adhesion on host cellular signaling and the outcome of bacteria-host interactions.

Recently a method was devised, whereby recombinant purified adhesins or adhesin fragments are directionally and covalently coupled to polymer particles similar in size to bacteria, thus mimicking bacterial surface display. This approach has been used on a range of different adhesins, from Gram-negative Multivalent Adhesion Molecules (MAMs)<sup>2, 3</sup>, to Gram-positive

adhesins including *Staphylococcus aureus* fibronectin binding protein (FnBPA)<sup>4, 5</sup> and *Streptococcus pyogenes* F1<sup>6, 7</sup>. This method has allowed the dissection of adhesin fragments important for host cell binding, identification of host surface receptors and determination of their behavior on the host cell surface, while taking both binding affinity and avidity into consideration<sup>1, 8</sup>. Additionally, this approach has been used to investigate the efficacy of immobilized adhesins and derivatives as inhibitors of host-pathogen interactions<sup>8, 9</sup>. It has been demonstrated that surface coupled derivatives of the *Vibrio parahaemolyticus* Multivalent Adhesion Molecule (MAM) 7 can be used to attenuate a range of bacterial infections *in vitro*, including those caused by multidrug-resistant pathogens such as *Acinetobacter baumannii* and methicillin-resistant *S. aureus* (MRSA)<sup>7, 9</sup>.

Herein, we describe different chemistries which can be used to immobilize adhesins to commercially available functionalized polystyrene particles. Due to the wide array of available surface functionalities, labels and particle sizes, these are a useful scaffold for the production of bacteriomimetic materials to investigate adhesin-host interactions. We further describe methods for the initial characterization of the coupling reaction, and for the calculation of important properties of the resulting materials. Finally, the use of bead-coupled adhesins as competitive inhibitors of *in vitro* bacterial infections is discussed as an example of their application, as well as the future scope and limitations of this approach.

## **PROTOCOL:**

### **1. Chemical coupling of proteins to polymer beads**

#### **1.1. Thiol-amine directional coupling**

Note: This protocol is suitable for coupling of cysteine containing proteins to amine-functionalized polymer beads, using Sulfosuccinimidyl 4-(*p*-maleimidophenyl)butyrate (Sulfo-SMPB) as cross-linking agent (Figure 1).

**[Place Figure 1 here]**

##### **1.1.1 Preparation of reagents:**

1.1.1.1 Prepare PBS (100 mM sodium phosphate, 150 mM NaCl, pH 7.0) and autoclave.

1.1.1.2 Prepare a 100x stock (0.5 M or 287 mg/ml in PBS) of TCEP (tris(2-carboxyethyl)phosphine) immediately before use.

1.1.1.3 Prepare a 5x stock (10 mM or 4.58 mg/ml in dH<sub>2</sub>O) of Sulfo-SMPB immediately prior to use.

1.1.1.4 Prepare a 10x stock (500 mM or 88 mg/ml in PBS, pH 7.0) of cysteine immediately before use.

133

134 **1.1.2 Bead activation:**

135

136 1.1.2.1 Mix the bead suspension by gently inverting and transfer the required amount of bead

137 suspension (e.g., 12  $\mu$ l) into a sterile 1.5 ml tube containing 1 ml sterile PBS, pH 7.0.

138 1.1.2.2 Gently pipette up and down to wash the beads and pellet by centrifugation in a

139 microcentrifuge (2 minutes at 16000 x g).

140 1.1.2.3 Carefully remove the supernatant with a pipette and discard.

141 1.1.2.4 Resuspend the bead pellet in 1 ml of fresh sterile PBS and repeat the washing step.

142 1.1.2.5 Resuspend the bead pellet in 0.8 ml of PBS.

143 1.1.2.6 Add 200  $\mu$ l of freshly prepared 10 mM Sulfo-SMPB, to give a final concentration of 2

144 mM.

145 1.1.2.7 Incubate the bead suspension for 1 hour at 25 °C on a rotating wheel.

146 **1.1.3 Protein reduction:**

147 1.1.3.1 During the incubation period of the activation step, prepare the protein for the

148 following coupling step so it can be immediately added to the activated beads.

149 Note: Although this reduction step is not always required for GST (Glutathione S-transferase)

150 fusion proteins, it is recommended to ensure a high coupling efficiency.

151 1.1.3.1.1 Check the protein concentration, and adjust it to the final concentration required

152 for the coupling reaction. Note: 6  $\mu$ M protein in PBS, and a volume of 1 ml are usually used.

153 1.1.3.1.2 Add TCEP stock to give a final concentration of 5 mM.

154 1.1.3.1.3 Incubate the solution for 30 minutes at room temperature.

155 Note: The reaction mixture can be directly used for the following coupling reaction.

156 1.1.3.1.4 Retain a small amount (a few  $\mu$ l) of the protein solution for determining the

157 protein concentration and calculation of coupling efficiency (see section 2).

158 **1.1.4 Protein coupling step:**

159 1.1.4.1 Pellet the activated beads by centrifugation (2 minutes, 16000 x g in a microcentrifuge),

160 and wash the pellet once in 1 ml of fresh sterile PBS.

161 1.1.4.2 Resuspend the pellet in the prepared protein solution (e.g., 1 ml), to give the desired

162 protein concentration.

Note: The protein concentration during the coupling step will depend on the average coupling efficiency (see section 2 for determination of the coupling efficiency) and desired coupling density (as calculated in 1.1.4.3). The efficiency is approximately 85% and the desired final concentration 5  $\mu$ M in the 10x bead suspension, so protein concentration during the coupling step should be approximately 6  $\mu$ M.

1.1.4.3 Calculate the coupling density using the following formula:

$$\text{coupling density } \rho_c = \frac{\text{protein conc}}{\text{protein } M_w} \cdot \frac{6 \cdot 10^{23}}{\text{bead conc}} \cdot \frac{1}{\pi d^2}$$

where

$\rho_c$  coupling density [number of protein molecules/nm<sup>2</sup>],  
 protein conc protein concentration [mg/ml],  
 protein  $M_w$  protein molecular weight [Da],  
 bead conc bead concentration [number of beads/L],  
 d bead diameter [nm]  
 $6 \cdot 10^{23}$  Avogadro's number

1.1.4.4 Use the coupling density to calculate the average ligand spacing:

$$\text{average ligand spacing [nm]} = \sqrt{\frac{1}{\rho_c}}$$

1.1.4.5 Incubate the protein-bead suspension for 2 hours at 25 °C on a rotating wheel.

Note: Some proteins may not be stable at room temperature. In these cases, the reaction can be carried out at 4 °C overnight.

1.1.4.6 Deactivate remaining activated groups on the beads by adding cysteine stock to a final concentration of 50 mM and incubate the suspension for 30 minutes at 25 °C on a rotating wheel.

1.1.4.7 Pellet the beads by centrifugation (2 minutes, 16000 x g in a microcentrifuge).

1.1.4.8 Keep the supernatant for determining the protein concentration and calculation of coupling efficiency (see section 2).

1.1.4.9 Wash the bead pellet twice with 1 ml PBS and resuspend in 1 ml fresh PBS to give the final product.

Note: The above protocol will typically give 1 ml of coupled protein, at a final concentration of 5  $\mu$ M protein, which can be used as a 10x stock for subsequent experiments (see section 3).

1.1.4.10 To proceed to section 3 of the protocol, work with 100  $\mu$ l/ml of bead stock, or at a final

concentration of 500 nM protein. Note: A good starting point for this procedure will be  $2 \cdot 10^{12}$  beads/ml, resulting in an average coupling density of  $3 \cdot 10^{-4}$  proteins/nm<sup>2</sup> or an average spacing of 57 nm on a bead of 2 µm diameter.

## **1.2. Thiol-carboxy directional coupling**

Note: This protocol is suitable for coupling cysteine containing proteins to carboxyl-functionalized polystyrene beads. The carboxyl moiety is first activated using 1-Ethyl-3-(3-dimethylaminopropyl)carbodiimide (EDC)/ N-hydroxysuccinimide (NHS), amine modified and then cross-linked using Sulfo-SMPB (Figure 2).

**[Place Figure 2 here]**

### **1.2.1 Preparation of reagents:**

1.2.1.1 Prepare PBS (100 mM sodium phosphate, 150 mM NaCl, pH 7.0) and autoclave.

1.2.1.2 Prepare a 100x stock (0.5 M or 287 mg/ml in PBS,) of TCEP immediately before use.

1.2.1.3 Prepare a 10x stock (20 mM, or 4 mg/ml in PBS) of EDC (1-Ethyl-3-(3-dimethylaminopropyl)carbodiimide) immediately before use.

1.2.1.4 Prepare a 10x stock (50 mM, or 6 mg/ml in PBS) of NHS (N-hydroxysuccinimide) immediately prior to use.

1.2.1.5 Prepare a 5x stock (10 mM or 4.58 mg/ml in dH<sub>2</sub>O) of Sulfo-SMPB immediately prior to use.

1.2.1.6 Prepare a 10x stock (500 mM or 88 mg/ml in PBS, pH 7.0) of cysteine immediately before use.

### **1.2.2 Bead activation:**

1.2.2.1 Mix bead suspension by gently inverting and transfer the required amount of bead suspension (e.g., 12 µl) into a sterile 1.5 ml tube containing 1 ml sterile PBS, pH 7.0.

1.2.2.2 Gently pipette up and down to wash the beads and pellet by centrifugation in a microcentrifuge (2 minutes at 16000 x g).

1.2.2.3 Carefully remove the supernatant with a pipette and discard.

1.2.2.4 Resuspend the bead pellet in 1 ml of fresh sterile PBS and repeat the washing step.

1.2.2.5 Resuspend the bead pellet in 0.8 ml of PBS.

235 1.2.2.6 Add 100 µl of 10x EDC stock (2 mM final concentration) and immediately 100 µl of 10x  
236 NHS stock solution (5 mM final concentration) to the bead suspension.

237 1.2.2.7 Incubate the bead suspension for 30 minutes at 25 °C on a rotating wheel.

238 1.2.2.8 Wash beads once in 1 ml of PBS and resuspend in 0.8 ml fresh sterile PBS.

239 1.2.2.9 Add 200 µl ethylenediamine, and incubate the bead suspension for 1 hour at 25 °C on a  
240 rotating wheel.

241 1.2.2.10 Wash beads once with 1 ml PBS and resuspend the beads in 0.8 ml of fresh PBS.

242 1.2.2.11 Proceed from section 1.1.2.6 (bead activation with Sulfo-SMPB) as described  
243 under section 1.1.2 and follow the rest of the protocol described in section 1.1 (Thiol-amine  
244 directional coupling). Perform protein preparation and protein coupling steps in a manner  
245 identical to those described in section 1.1.

246 Note: The typical coupling efficiency using this protocol is slightly lower (approx. 75%)  
247 compared to section 1.1, therefore adjust the initial protein concentration accordingly to  
248 achieve the same coupling density.

## 249 **2. Determination of coupling efficiency**

250 Note: To determine protein concentration, use Bradford Reagent <sup>10</sup> and colorimetric assays as  
251 follows:

252 2.1 Gently invert Bradford Reagent to ensure homogeneity of the reagent.

253 2.2 Using a 10 mg/ml BSA stock solution, prepare protein standards covering concentrations  
254 from 0.1 to 1.5 mg/ml BSA in buffer.

255 2.3 Add 250 µl of Bradford Reagent to wells of a 96-well plate. Prepare enough wells for all  
256 samples, protein standards and negative controls (buffer only).

257 2.4 Add 5 µl of protein sample (protein standards or buffer, for the negative control) to the  
258 reagent in the 96-well plate.

259 2.5 Incubate the plate on an orbital shaker at room temperature for 10 minutes.

260 2.6 Measure absorbance at 595 nm using a plate reader.

261 2.7 Generate a standard curve of BSA concentration versus A<sub>595nm</sub> and use this to  
262 determine protein concentrations in initial and supernatant samples.

263 2.8 Calculate the concentration of coupled protein as follows:  
264



$$[Coupled\ protein] = [Initial\ protein] - [Supernatant\ protein]$$

2.9 Calculate the coupling efficiency as:

$$Coupling\ efficiency\ (\%) = \frac{[Initial\ protein] - [Supernatant\ protein]}{[Initial\ protein]} \times 100$$

### 3. Use of bead-coupled adhesins in competition assays

#### 3.1 Preparation:

3.1.1 Seed 1 ml per well of Hela cells at a concentration of 150000 cells/ml into a 24-well plate the day before the competition assay, to allow cells to reach approximately 80% confluency prior to starting the experiment.

3.1.2 Set up each experimental condition in triplicate. Include wells for negative controls (no bacteria added during competition assays), positive controls (control beads coupled to fusion-tag only added during competition assays) and lysis controls (for cytotoxicity experiments).

3.1.3 Inoculate a 5 ml marine LB (MLB) culture with a fresh colony of *V. parahaemolyticus* and grow overnight at 30 °C, shaking.

3.1.4 Prepare sufficient bead-coupled MAM, as described in section 1. (Coupling). Allow for 100 µl of 10x bead stock per well.

#### 3.2 Competition assay:

3.2.1 On the day of the competition experiment, measure the OD<sub>600</sub> of the bacterial culture.

3.2.2 Prepare infection media by diluting bacterial cultures into colorless DMEM without additives, pre-warmed to 37 °C, containing 10% v/v bead suspension (either adhesin-coupled beads or control beads), to give an MOI of 10. Prepare 1 ml/well and 10-20% excess volume per sample.

Note: For the above mentioned conditions (24-well plate, *V. parahaemolyticus*, MOI 10), the necessary volume of overnight culture to be added per well (µl/ml) is calculated as 3/OD<sub>600</sub> of the culture. For example, 1 ml of infection medium will typically contain 100 µl bead suspension, a few µl of bacterial culture and be made up to 1 ml with colorless, unsupplemented DMEM.

3.2.3 Remove old medium from wells and wash cultured Hela cells by adding 1 ml of sterile PBS pre-warmed to 37 °C to each well.

3.2.4 Remove PBS and add 1 ml of infection medium per well. Also set up controls, by adding solutions containing control beads and bacteria (positive control) or adhesin beads and no

298 bacteria (negative controls), or DMEM containing 0.1% Triton X-100 (lysis control, only  
299 necessary for cytotoxicity measurements).

300 3.2.5 Incubate the plate in a tissue culture incubator at 37 °C for the desired amount of time  
301 (e.g. 4 hours for cytotoxicity measurements or 1 hour for adhesion measurements).

302 Note: Both bacterial adhesion and cytotoxicity on host cells can be used as read-outs for the  
303 efficacy of inhibition. If the time points of cytotoxicity and adhesion measurements coincide,  
304 both assays may be performed using samples from the same well, since cytotoxicity is  
305 determined using the culture supernatant and attachment assays use samples derived from the  
306 remaining cell layer.

### 307 **3.3 Cytotoxicity measurements:**

308 3.3.1 At indicated time points (e.g., 4 hours post infection), remove three times 200 µl from  
309 each 24-well and transfer to a 96-well plate.

310 3.3.2 Spin 96-well plates at 1500 x g, 5minutes and transfer 100 µl from each well into a fresh  
311 96-well plate.

312  
313 3.3.3 Add 100 µl of the media used during infection experiments to fresh wells of the 96-well  
314 plate in triplicate (these will be used as blanks).

315  
316 3.3.4 Carry out the lactate dehydrogenase (LDH) release assay using a LDH cytotoxicity  
317 detection kit and following the manufacturer's instructions.

318  
319 3.3.4.1 Briefly, calculate the amount of reagent needed in increments of 25 (e.g., if 62 samples  
320 are to be measured, make up enough reagent mix for 75 etc.).

321  
322 3.3.4.2 For example, for 100 samples, mix 11.25 ml of reagent A with 250 µl of reagent B.  
323 Invert, do not vortex to avoid foaming.

324  
325 3.3.4.3 Put the mixture in a reservoir to be able to pipet with a multi-channel pipette.

326  
327 3.3.4.4 Add 100 µl of reagent mix to each sample.

328  
329 3.3.4.5 Incubate plate at room temperature and read the absorbance at 490 nm on a plate  
330 reader at 10, 20, 30 minutes.

331  
332 3.3.4.6 Analyze the data set for which the absorbance of the lysis control sample is high but still  
333 within the linear range of the plate reader (typically, 2-3 absorbance units).

334  
335 3.3.4.7 Express results as % cytotoxicity, using the following formula for conversion:

336

$$\text{Cytotoxicity (\%)} = \frac{A490 (\text{sample}) - A490 (\text{blank})}{A490 (\text{triton control}) - A490 (\text{blank})} \cdot 100$$

### 3.4 Measurement of bacterial adhesion:

3.4.1 At indicated time points (e.g., 1 hour post infection), remove media from the cell layer.

3.4.2 Thoroughly wash the cell layer with sterile, pre-warmed PBS (at least 3-4 washes of 1 ml PBS each) to remove any un-attached cells.

3.4.3 Lyse host cells by adding 1 ml of a sterile 1 % v/v Triton X-100 solution in PBS per well.

3.4.4 Incubate the plate at 37 °C for 5 minutes.

3.4.5 Pipette each sample up and down several times before transferring the contents of each well to separate 1.5 ml tubes.

3.4.6 Prepare 10-fold serial dilutions of samples into sterile PBS (e.g., use 100 µl of sample and 900 µl of PBS).

3.4.7 Plate 100 µl of each sample on MLB agar and spread using a cell spreader. Optimize which dilutions to plate depending on the bacterial strains and time point. Note: For the described experimental setup, 10<sup>5</sup> or 10<sup>6</sup> fold dilutions give a suitable number of CFUs.

3.4.8 Incubate plates at 37 °C overnight and enumerate bacteria by colony counting.

### REPRESENTATIVE RESULTS:

The *V. parahaemolyticus* adhesin MAM7 contains seven tandem mammalian cell entry (MCE) domains involved in recognition of host surface receptors. We used polystyrene beads coupled to recombinant, purified fragments encompassing either all seven tandem MCE domains (MAM7) or only the first MCE domain (MAM1) to test the ability of these materials to compete with *V. parahaemolyticus* for host cell binding and the resulting efficacy of these materials as adhesion inhibitors. Hela cells were infected with *V. parahaemolyticus* strain POR1, and cytotoxicity resulting from *in vitro* infection was evaluated after 4 hours (Figure 3). Treatment of Hela cells with 0.1% Triton X-100 (positive lysis control) resulted in complete cell lysis, uninfected cells displayed very low levels of cytotoxicity. *In vitro* infection of untreated cells with POR1 resulted in very high levels of cell lysis, and this was inhibited by MAM7-coupled beads but not MAM1-coupled or GST-coupled control beads (Figure 3).

### [Place Figure 3 here]

Enumeration of *V. parahaemolyticus* attached to either untreated Hela cells, or cells incubated with MAM7-, MAM1-, or GST control beads, revealed that MAM7-beads but not MAM1- or GST- control beads outcompete *V. parahaemolyticus* for attachment to host cell surface


receptors (Figure 4).

[Place Figure 4 here]


Adhesins coupled to fluorophore labeled beads result in the generation of materials that mimic bacterial adhesion to host cells and are powerful tools for cellular imaging (Figure 5). Using MAM7-coupled fluorescent blue beads, we characterized the process of MAM7-mediated attachment to epithelial cells. Attachment of MAM7 to host cells resulted in actin rearrangements and formation of stress fibers, which were co-visualized using rhodamine-phalloidin to stain for F-actin (Figure 5B).

[Place Figure 5 here]


### Figures & Legends:




**Figure 1: Cross-linking strategy used for directional thiol-amine coupling of proteins to polymer beads.** Amine-modified polystyrene beads are activated with Sulfo-SMPB. The maleimide reacts with free cysteines to directionally couple proteins to beads.




**Figure 2: Cross-linking strategy used for directional thiol-carboxyl coupling of proteins to polymer beads.** Carboxylated polystyrene beads are first activated with EDC and modified with NHS to form a semi-stable NHS-ester. Subsequent amine-coupling of ethylenediamine results in a free amine group, which reacts with Sulfo-SMPB. Maleimide activated beads can then react with free cysteines to directionally couple proteins to beads.



**Figure 3. Characterization of *Vibrio parahaemolyticus* induced cytotoxicity in epithelial cells during competition assays.** Cells were treated with *V. parahaemolyticus* (V.p), indicated with (+), or no bacteria (-) in the presence of different competing entities (comp.) as indicated. These included either no comp. (-), MAM7 beads (MAM7), MAM1 beads (MAM1) or GST control beads (GST). As controls, cells were treated either with Triton X-100 (triton, lysis control), or left uninfected (uninf). LDH release after 4 hours was measured as described in section 3, and results normalized to Triton controls (100 %) and blanks (0 %) as described above. Results are means  $\pm$  s.e.m from triplicate experiments.



**Figure 4. Characterization of bacterial attachment during competition experiments.** Hela cells were infected with *V. parahaemolyticus* POR2 (V.p. +), in the absence (-) or presence of competing entities (comp.) as follows: MAM7 beads (MAM7), MAM1 beads (MAM1) or GST control beads (GST). Bacterial adhesion was measured after 1 hour, as described in section 3. Results are means  $\pm$  s.e.m from triplicate experiments. Means (f.l.t.r. in CFU/ml) are  $2.30 \cdot 10^6$ ,  $1.53 \cdot 10^5$ ,  $2.05 \cdot 10^6$ ,  $2.12 \cdot 10^6$ .



**Figure 5. Preparation and use of fluorophore labeled biomimetic beads for imaging purposes.** (A) Suspension of fluorescent blue biomimetic beads (right tube, 10x stock) and buffer control

(left tube). **(B)** Attachment of MAM7-coupled beads to Hela cells results in actin rearrangements and stress fiber formation. Attachment of fluorescent blue beads and resulting actin stress fibers (red, stained with rhodamine-phalloidin) were imaged by microscopy. Bar, 10  $\mu$ m. Images in panel B were adapted from Lim *et al.*<sup>1</sup> and reproduced under the Creative Commons Attribution license.

## DISCUSSION:

Herein, we describe two protocols, which can be used to couple thiol-containing proteins to amine- or carboxylate-modified polystyrene beads, respectively. Due to ease of the procedure, thiol-amine coupling is preferable, but depending on the desired bead specification (diameter, fluorescence properties), use of amine-functionalized beads may not be possible and we have therefore included a protocol which will convert the carboxylate- into an amine moiety to give the researcher the greatest possible flexibility in choice of scaffold. Although both thiol-amine and thiol-carboxyl coupling work for any cysteine containing protein, directional coupling (i.e., immobilization of the protein per its N-terminus, mimicking surface display) requires a protein without cysteine residues, that is produced as a GST fusion protein, or that contains a single terminal cysteine residue introduced by site directed mutagenesis. If multiple reactive cysteines are contained within the protein, this will lead to random immobilization which may impede protein function. Many bacterial adhesins do not contain cysteines naturally. For others, these may be removed by site directed mutagenesis, although this would require extensive assays to ensure native structure and function are retained in the mutant. For GST fusion proteins, purified GST-tag coupled to beads can be used as a suitable negative control. Using uncoupled beads as a control should be avoided, as these often have a higher tendency to clump together or adhere to cells non-specifically. A simplified version of protocol 1.2., using only EDC, can be used to couple proteins to carboxyl-functionalized polymer beads, however in this case coupling takes place via primary amines in the protein and therefore does not guarantee directional coupling.

TCEP as a reducing agent must not be replaced with other commercially available and commonly used reducing agents, such as dithiothreitol (DTT) or 2-mercaptoethanol (BME), as the thiols contained within them will compete with protein coupling in the thiol-maleimide coupling step. PBS may be replaced with other buffers, but with the following considerations: Buffers may not contain primary amines (so Tris-containing buffers are not suitable). Use of buffers containing very low (< 10mM) salt concentrations leads to bead clumping and should also be avoided. Protein purity is also an important factor to consider during this procedure, and to achieve high quality data, pure proteins should be used. We routinely purify proteins in multiple steps, including at least an affinity purification and gel filtration step, but in some cases ion exchange chromatography is done as a third step. As a result the purity of proteins used for coupling is usually 90% or higher, as judged by SDS-PAGE.

It is recommended to determine protein concentrations both in the initial reaction mixture as well as of the supernatant after reaction completion. This will help to determine the apparent concentration of bead-coupled protein, and thus coupling density. Determination of both

values will also allow calculation of the coupling efficiency. This can be taken into account when preparing the initial protein solution in subsequent reactions, to achieve the desired final concentration and coupling density. Bradford reagent is particularly suited for determination of protein concentration before and after the coupling reaction, as none of the substances in the reaction interferes with the dye complex formation at the concentrations used. If lower protein concentrations are to be used, this method may have to be replaced by a more sensitive detection method, however attention has to be paid to the fact it has to be compatible with the substances contained within the coupling reaction. It is also recommended to use freshly prepared reagents and handle powdered reagents with care (eg, store in a sealed container and use silica beads, to avoid the reagents drawing moisture) since the quality of reagents will influence the coupling efficiency. If the coupling efficiency is lower than expected, possible remedies include increasing the initial bead and protein concentrations. If higher concentrations are being used, the concentration of coupling reagents has to be increased proportionally to ensure sufficient molar excess. Modifying bead/protein concentrations is usually a better step towards optimization rather than increasing reaction times. Since the protocols for bead coupling are lengthy, we commonly prepare a large batch of material. Aliquots of the suspension can be snap-frozen in liquid nitrogen and stored at -20 °C for several months. Thawed aliquots should not be refrozen and should be kept at 4 °C and used within 1-2 days. However, this will vary with the nature and stability of the protein used as should be tested on a small batch initially.

Bead-coupled adhesins can be used for many applications, as discussed below. This protocol described an assay that is commonly used to measure inhibition of bacterial binding and pathogen-mediated cytotoxicity on host cells. The assay is commonly performed to measure the capacity of bead-coupled MAMs to competitively inhibit infection of Hela epithelial cells with the sea-food borne pathogen *Vibrio parahaemolyticus*, using either a decrease in bacterial attachment to host cells or reduced cytotoxicity as a read-out. In both cases, preparation and competition assays follow the same protocol. Depending on the readout, different strains of *V. parahaemolyticus* are being used: the cytotoxic strain POR1 is used for cytotoxicity measurements, while the non-cytotoxic strain POR2 is used for measuring bacterial adhesion, since cell death and cell detachment compromises the procedure for quantifying attached bacteria.

Initially, competition experiments were set up as a step-wise protocol, where host cells were first pre-incubated with beads prior to the addition of bacteria. For *V. parahaemolyticus* and the bead specifications used (2 µm beads coupled to MAM7), both beads and bacteria can be added at the same time without changes in the resulting cytotoxicity. I.e., in this experimental setup, beads outcompete bacteria for host cell binding. Depending on the bacterial species and bead geometry used, there may be good reasons for maintaining the bead adhesion and bacterial infection as two separate steps. For example, to infect cells with non-motile bacteria, plates are commonly centrifuged after addition of the infection media. However, centrifugation of plates containing bead suspensions should be avoided, since this leads to a highly uneven distribution of beads on the cell layer. If smaller particle sizes are being used, beads will take longer to settle on the cell surface, in which case sufficient time should be allowed for bead

attachment prior to the infection. When bacterial adhesion is used as a read out, samples should be taken at time points where host cells are not significantly damaged by the infecting strain, as cell detachment and lysis can compromise the quantification of attached bacteria.

Instead of enumerating bacterial adhesion by dilution plating, samples may alternatively be processed for imaging (Figure 5). In this case, tissue culture cells should be seeded onto glass cover slips, rather than directly into wells. Additionally, fluorescent beads and bacteria expressing a fluorescent protein may be used, along with infection-specific host cell markers. For example, competition experiments are commonly imaged using fluorescent red rhodamine-phalloidin to stain the host cells' actin cytoskeleton and assess morphological changes resulting from infection, together with fluorescent blue beads and fluorescent green (GFP expressing or SYTO18 stained) bacteria.

A range of bead-coupled adhesins, including *Staphylococcus aureus* FnBPA, *Streptococcus pyogenes* F1 FUD and *Vibrio parahaemolyticus* MAM, have been used as biomimetic materials to study adhesion, adhesion inhibition and the contribution of adhesion to pathogen-mediated cytotoxicity<sup>1, 7, 8</sup>. One of the advantages of using this approach is the ease of visualization of attachment events, since the polymer beads used as scaffolds are available in a wide range of colors (e.g., blue, fluorescent red, blue, green, orange). Thus, direct protein labeling, which may interfere with function, can be avoided. Additionally, surface coupling mimics the multivalent display of adhesins on the bacterial surface, thus reflecting a more physiologically relevant conformation compared to soluble proteins.

Compared to studies using intact bacteria or bacterial mutants, the bead approach circumvents problems associated with bacterial growth. For example, longer-term (e.g., overnight) studies of bacterial adhesion to host cells using intact bacteria are often compromised by phenomena accompanying bacterial growth – acidification of the growth medium and nutrient depletion negatively affect host cells, and bacterial replication eventually compromises imaging quality.

More recently, the use of bead-coupled adhesins has been extended to include their use as tools for affinity purifications of host cellular factors involved in signaling processes downstream of bacterial attachment. *V. parahaemolyticus* MAM7, via binding to phosphatidic acids in the host cell membrane, triggers RhoA activation and actin rearrangements<sup>1, 3</sup>. MAM-coupled beads are being used to purify and identify proteins involved in the signaling platforms assembled as a consequence of MAM-host cell binding. Since beads can easily be separated from the supernatant by a short centrifugation step, and the protein of interest is covalently coupled, this is a good method to achieve separation from contaminant proteins and enrich relevant protein complexes, which can be used for downstream applications such as proteomics or Western Blotting.

#### ACKNOWLEDGMENTS:

The authors would like to thank members of the Krachler group for critical reading of the manuscript. DHS, DV and AMK were funded by the Biotechnology and Biological Sciences



Research Council (BB/L007916/1), FA was funded by a Republic of Iraq Ministry of Higher Education and Scientific Research Scholarship and NPS was funded by a CONICYT Scholarship.

## DISCLOSURES:

The authors have nothing to disclose.

## REFERENCES

1. Lim, J., Stones, D. H., Hawley, C. A., Watson, C. A. & Krachler, A. M. Multivalent Adhesion Molecule 7 clusters act as signaling platform for host cellular GTPase activation and facilitate epithelial barrier dysfunction. *PLoS Pathog* **10** (9), e1004421, doi:10.1371/journal.ppat.1004421 (2014).
2. Krachler, A. M., Ham, H. & Orth, K. Outer membrane adhesion factor multivalent adhesion molecule 7 initiates host cell binding during infection by gram-negative pathogens. *Proc Natl Acad Sci U S A* **108** (28), 11614-11619, doi:10.1073/pnas.1102360108 (2011).
3. Krachler, A. M. & Orth, K. Functional characterization of the interaction between bacterial adhesin multivalent adhesion molecule 7 (MAM7) protein and its host cell ligands. *J Biol Chem* **286** (45), 38939-38947, doi:10.1074/jbc.M111.291377 (2011).
4. Joh, D., Speziale, P., Gurusiddappa, S., Manor, J. & Hook, M. Multiple specificities of the staphylococcal and streptococcal fibronectin-binding microbial surface components recognizing adhesive matrix molecules. *Eur J Biochem* **258** (2), 897-905, doi:10.1046/j.1432-1327.1998.2580897.x (1998).
5. Bingham, R. J. *et al.* Crystal structures of fibronectin-binding sites from *Staphylococcus aureus* FnBPA in complex with fibronectin domains. *Proc Natl Acad Sci U S A* **105** (34), 12254-12258, doi:10.1073/pnas.0803556105 (2008).
6. Ensenberger, M. G. *et al.* Specific interactions between F1 adhesin of *Streptococcus pyogenes* and N-terminal modules of fibronectin. *J Biol Chem* **276** (38), 35606-35613, doi:10.1074/jbc.M105417200 (2001).
7. Hawley, C. A., Watson, C. A., Orth, K. & Krachler, A. M. A MAM7 peptide-based inhibitor of *Staphylococcus aureus* adhesion does not interfere with in vitro host cell function. *PLoS One* **8** (11), e81216, doi:10.1371/journal.pone.0081216 (2013).
8. Krachler, A. M., Ham, H. & Orth, K. Turnabout is fair play: use of the bacterial Multivalent Adhesion Molecule 7 as an antimicrobial agent. *Virulence* **3** (1), 68-71, doi:10.4161/viru.3.1.18172 (2012).
9. Krachler, A. M., Mende, K., Murray, C. & Orth, K. In vitro characterization of multivalent adhesion molecule 7-based inhibition of multidrug-resistant bacteria isolated from wounded military personnel. *Virulence* **3** (4), 389-399, doi:10.4161/viru.20816 (2012).
10. Bradford, M. M. A rapid and sensitive method for the quantitation of microgram quantities of protein utilizing the principle of protein-dye binding. *Anal Biochem* **72**, 248-254, doi:10.1016/0003-2697(76)90527-3 (1976).